

Case Study

JULABO ED-39 and FT402

Determining the lowest working temperature when using an FT402 immersion cooler with an ED-39 heating bath circulator

Case Study
03-2007
Page 1 from 1



The Objective

The test is carried out to determine the lowest temperature that could be achieved when using an FT402 for cooling an ED-39

Setup


For the test an ED-39 heating bath circulator (filling volume: 39 liters) is filled with Thermal G bath fluid. The FT402 immersion cooler is used as cooling device; the probe is immersed at the left side of the bridge. The temperature inside the bath is registered with an LC6 laboratory controller connected to a laptop.

Test Conditions

| | |
|---------------------|---|
| Instruments | ED-39 heating bath circulator FT402 immersion cooler |
| Mains voltage | 230 Volts/ 50 Hz |
| Ambient temperature | 20°C (room temperature) |
| Bath fluid | Thermal G |

Result


Under the test conditions above, a value of -9.9 °C is achieved as lowest temperature.



39 liters bath with immersion probe for cooling the bath



Setup with heating bath circulator and immersion cooler



Cool-down curve

JULABO GmbH
Eisenbahnstraße 45 • 77960 Seelbach / Germany
Tel. +49 (0) 7823 51-0 • Fax +49 (0) 7823 24 91
info@julabo.de • www.julabo.de