

## CORIO C-BT19 Open Heating Bath Circulator

**The new generation of laboratory circulators for routine work and standard tasks**

The Open Heating Bath Circulators of the new CORIO™ series feature high-quality transparent bath tanks.

### Your advantages

- Models for internal applications
- Bright, white, easy to read display
- Very quiet



**Made**  
in Germany

### Technical Data

Order No.	9011319
Model series	CORIO
Category	Open Heating Bath Circulators
Working temperature range (°C)	+20 ... +100
Temperature stability (°C)	±0.03
Setting / display resolution	0.01 – 0.1 °C
Temperature Display	LED
Heating capacity (kW)	1
Circulation capacity Flow rate (l/min)	6
Circulation capacity pressure (psi)	1.45
Bath opening / bath depth (W x L / D inch)	11.8 x 13.7 / 5.9
Filling volume liters	14 ... 19
Ambient temperature	5...40 °C
Dimensions W x L x H (inch)	15 x 22.9 x 15
Weight (LBS)	15.5
Classification according to DIN12876-1	Classification I (NFL)
Cooling coil	optional
Bath tank	Polycarbonat
Available voltage versions	208 - 230 V / 50 - 60 Hz 100 - 115 V / 50 - 60 Hz

## Characteristics



### 100% Checked.

100% testing. 100% quality. Each CORIO™ undergoes thorough quality testing before leaving the factory.



### Modern. Reliable.

High-grade components in every CORIO – platinum sensors, proven motor technology, CAN-Bus communication and much more.



### Brilliant.

Bright display with strong illumination can be viewed easily from a long distance.



### Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



### Unique.

Specially designed pyramids ensure that the condensed liquid drips directly back into the bath.



### Safety.

CORIO CD complies with Class III (FL) according to DIN 12876-1 and switches off automatically in case of high temperature or low liquid level alarm.



### Green technology.

Development consistently applied environmentally friendly materials and technologies.



### Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



### JULABO. Quality.

Highest standards of quality for a long product life.



### Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at [www.julabo.com](http://www.julabo.com).



### Locked in.

The lockable power plug guarantees safe connection. More process safety.



### Switch on. And off you go.

Intelligent operating concept. Ready for operation with just a few quick and easy steps.